

FERNO 419 BABY CAPSULE RESTRAINT

User Instructions

The 419 Baby Capsule Restraint is an Australian Standards approved harness that provides safe restraint for rearward facing baby capsules on stretchers.

USER INSTRUCTIONS

1. Remove any existing harness or restraint straps so the stretcher mattress is clear with no obstructions for fitting baby capsule and restraint




2. Fit the baby capsule outer shell on the stretcher and raise the backrest so it aligns with the baby capsule outer shell.




3. Attach the smaller strap to the stretcher main frame harness mount point. Extend the webbing to full length in preparation for fitting



4. Feed the unattached end of the smaller strap through the underside of the baby capsule outer shell and thread through to the opposite side.



5. Attach to the harness main frame mount point on stretcher and pull excess webbing to tighten strap to baby capsule outer shell.



6. Place the baby capsule inside the outer shell.

7. Fit the two large metal clips from the large restraint to the moulded sides at the rear of the outer shell and push up to lock tight.



FERNO 419 BABY CAPSULE RESTRAINT User Instructions

8. Run the quick release buckles to the stretcher backrest harness mount points and pull excess webbing tight to adjust the capsule into the stretcher backrest.



9. Take the baby capsule tether strap and loop around the foot end of the stretcher main frame.



10. Attach the baby capsule tether strap to the standard baby capsule strap and tighten.



11. Check the baby capsule is firmly restrained and all straps are tight. Apply baby capsule restraints to the baby as per standard capsule protocols.

The baby capsule is now ready to safely transport baby.



CLEANING & MAINTENANCE

- Yearly maintenance inspections recommended and should be complete by a competent person
- Alcohol/ bleach based products not recommended for cleaning restraints
- Restraint webbing can be immersed in approved cleaning solution. DO NOT immerse buckles in liquid. Dry in shade.
- Wipe/ spray metal components with approved cleaning products following manufacturer's instructions. DO NOT immerse buckles in liquid.
- Restraints which are removed from a stretcher should be visually inspected before re-fitting to the stretcher.

Examine restraints for:

- Cuts/ abrasions in webbing and frayed/ broken stitching
- Oil, acid or alkaline contamination
- Wear, distortion, fractures or corrosion in the metal fittings or buckles
- Correct functioning of buckles
- Correct fitment of straps in buckles.

CUSTOMER SERVICE: Toll free in Australia 1800 804 647
E: info@ferno.com.au W: www.ferno.com.au