



Ferno Patient Slides

DHS-485A / DHS-485 / DHS-486 / FWE-SPS



This range of traditional Ferno patient slides is available in three sizes, including bariatric, for lateral transfers of supine patients to and from hospital beds, stretchers and similar.

The black FWE-SPS Patient Slide is designed to QAS specifications.

Specifications	DHS-485A	DHS-485	DHS-486	FWE-SPS
Description	Small (pair)	Standard	Large	QAS Specification
Colour	White	White	White	Black
Length	750 mm	1650 mm	1611 mm	970 mm
Width	570 mm	495 mm	640 mm	457 mm
Depth	10 mm	10 mm	10 mm	2.5 mm
SWL	160 kg	160 kg	160 kg	
Material	HDPE	HDPE	HDPE	HDPE