

(ANZ VERSIONS) **POWERX** POWER STRETCHER

DATA SHEET



POWERX AU Multi-Purpose Stretcher (FWE-POWERX)

Extra Features Included

Integrated IV / Push Poles (2) plus

pole pads (2)

Monitor Hook behind headrest

Adjustable and detachable winged headrest

Shoulder padding to patient harness for added comfort

Extender Pillow for SX with Extender Wings

SX Surface Extender with Extender Wings

ABPR Harrier Mattress with hook and loop fastenings

Anti-static Wheels (4)

Multi-Directional Wheels (4) each with individual brake and footactivated directional wheel lock

Winch Attachment

Stow-Net

CD/ZA/MGS Cylinder Holder

Designed for use in ambulances with a tail-lift or ramp, POWERX is ergonomically designed for enhanced handling and manoeuvrability, making it easier and safer to operate.

Powered to minimise potential manual handling injuries, the feature-rich POWERX has been further configured to satisfy the needs of paramedics, bariatric and specialist caregivers in Australia and New Zealand.

EASY TO USE & MANOEUVRE

- Simple push button controls are intuitive to use and learn
- Four swivelling anti-static wheels each with an individual brake and foot-activated directional wheel lock
- Four push handles at each corner of the stretcher. Two of these push handles incorporate extendable IV pole tips with hook and loop fastenings for fluid bags
- The POWERX shortens to 1450mm or 1540mm with the Stow-Net installed
- Four bump wheels on each stretcher corner prevents potential damage to building surfaces while protecting the stretcher from accidental damage

PATIENT COMFORT & SAFETY

- The mattress is constructed from layers of high density and memory foam and covered in a robust 2-way stretch material coated on both sides to repel fluids.
- Flat mattress surface facilitates unhindered lateral transfers
- The winged headrest can be adjusted to cradle the patient's head for enhanced comfort
- An upgraded multi-point patient harness and leg restraint set includes customised shoulder padding for extra patient security and comfort
- New SX with Extender Wings and cushioned pads increase surface area

EASY CLEANING & MAINTENANCE

1

- The stretcher's moulded top bed surface is light in weight but durable and easy to clean and maintain
- The mattress and headrest are both removable and offers enhanced infection control with bacteria resistant, wipe-clean surfaces.

POWER CONTROL

- Concens Controller and 28V (5Ah)
 Milwaukee battery
- Battery delivers ±110 lift cycles from one charge based off a 200kg patient, and ±66 lift cycles at maximum SWL of 350kg.
- Short recharge time (<1hr), autosleep mode activates after 1 minute of inactivity





POWERX AU Standard Stretcher (POWERX)



POWERX AU Multi-Purpose Stretcher (FWE-POWERX)



POWERX AU Multi-Purpose Stretcher with optional PacRac+ Extension Kit and PacRac+

Stretcher Inclusions (AU Conversions)

DESCRIPTION	Standard	Multi-Purpose
	POWERX	FWE-POWERX
FWE-POWERX Multi-Purpose Stretcher with two-part Stainless Steel Lock		FWE-POWERXL
Monitor Hook	•	•
CD/ZA/MGS Cylinder Holder	•	•
Winch-Link Attachment	•	•
Headend Storage Net	•	•
Adjustable ABPR winged headrest	•	•
Standard 2-Piece 5' Patient Restraints	•	_
Blue Patient Harness & Leg Restraints (BRH-417-POWERX)	_	•
Padded Wide Shoulder Pads (BRH PADDING)	_	•
ABPR (Standard) Mattress with anti-bacterial VELCRO® Fasteners	•	_
ABPR (Harrier) Mattress with anti-bacterial VELCRO® Fasteners	_	•
Side LED Lighting Strips	•	•
Anti-Static Wheels	4	4
Directional Wheel Locks	4	4
Wheel Brakes	4	4
Standard Push Handles	4	2
Integrated IV / Push Handles (90-0364)	_	2
Push Handle Padding	_	2
SX Surface Extenders (FWESX)	2	2
SX Surface Extender Pillow Pads	2	_
SX Extender Wings (FWESXEW)	_	2 sets
SX Extender Pillow Pads (FWESXEP)	_	1 pair
AS/NZ: 4535 Upgrade Kit and Backrest	_	•
240V Mains Charger for Milwaukee Battery	1	1
12V Battery Charger for Milwaukee Battery	_	1
Milwaukee 28V Battery	2	2

Stretcher Options

PacRac+ Extension Kit for PowerX (009010150)

PacRac+ Equipment Table (FWE-PR+ ITX)

Specifications

DESCRIPTION	Standard	Multi-Purpose
LENGTH: Shortened / Extended Shortened with Stow-Net fitted	1450 / 1930 mm 1540 mm	
HEIGHT: Low / Raised	480 / 880 mm	
WIDTH: Normal With both SX horizontally positioned	570 mm 920 mm	
SWL: Raised / Lowered	350 kg / 400 kg	
POWER:	28 Volt DC Battery 5Ah	
WEIGHT:	≈88 kg	≈92.5 kg
FASTENER: (2-Part Lock)	N/A	13 kg
STRETCHER POSITIONING:	TrendelenburgShock PositionKnee Contour / Gatch	
LED SIDE LIGHTING OPTIONS:	Constant white / Flashing white / Flashing red / Alternate flashing red and white	