



## iNCASE Portable Storage



The new iNTRAXX iNCASE is a semi-rigid case designed with integrated pouch mount carriers to securely hold SafePouches inside the case. Strong yet lighter than its predecessor, the iNCASE is a streamlined upgrade from the original SafePak, which has now been retired.

Contents remain visible and easily accessible while the case is mounted on a vehicle wall, but can be quickly removed for use in the field. The iNCASE and SafePouch system is ideal for “make-ready” or swap-n-go restocking methods to efficiently resupply vehicles.

### Features

- Semi-rigid casing is softer and safer in the event of a head strike during an accident
- User-friendly rubberised handles lock together for easy carrying
- Heavy-duty zip with glove-friendly pulls
- Easy to clean and decontaminate
- Supplied as requested without pre-configured pouches
- Pouches and Pouch Kits are sold separately to suit your specific needs



Heavy-duty zip with glove-friendly pulls




Interlocking rubberised handles for easy carrying

<b>Product Part #</b>	0480252
<b>Overall Size: mm / inches</b>	(L) 455 x (W) 467 x (D) 252 mm (L) 17.9 x (W) 18.4 x (D) 9.9 in.
<b>Weight</b>	4.9 kg / 10.8 lbs
<b>Load Capacity</b>	12 kg / 26.5 lbs
<b>Rating</b>	10g / 22g / 26g



Configure iNCASE with your preferred SafePouch layout



### Glossary

W: Width

D: Depth

L: Length